

SAFETY

SHOP www.belart.com



Emergency Eye Wash Safety Stations

Built-In Eye Cup Provides Immediate Emergency Care

U.S. Patent 6,540,726

This fluid delivery system features a built-in eye cup that allows the injured person to be treated in any position, whether standing, sitting, or lying down and provides effective and immediate eye irrigation. Quick instruction diagrams on the bottle and station confirm how to use the eyewash bottle, and multi-language directions are on the station (English, French, German and Spanish). Three stations are available to suit all workstations: a single or a double station with 1000ml (32 oz) eyewash bottle(s) or a single station with a 500ml (16 oz) bottle. Stations are shipped with empty bottles and can be filled with F24879-0032 Sterile Eye Wash Solution Refills (page 311). Maintenance instructions for bottles are included. Single station dimensions: 44.5W x 38cmH (17½ x 15"). Double station dimensions: 45.7W x 34.9cmH (18 x 13¾"). **Per each, 12 per case.**

CATALOG NO.	DESCRIPTION
F24865-0000	Single Station with 500ml (16 oz.) Empty Eye Wash Bottle
F24866-0000	Single Station with 1000ml (32 oz.) Empty Eye Wash Bottle
F24868-0000	Double Station with two 1000ml (32oz) Empty Eye Wash Bottles
F24850-0000	500ml (16 oz.) Empty Eye Wash Bottle
F24851-0000	1000ml (32 oz.) Empty Eye Wash Bottle



Built-in Eye Cups



Works in Any Position



All Scienceware® and PLUM Emergency Eye Wash Products Serve as Immediate First Aid only. Further Examination and Treatment by a Physician are Required.